



Crile Hemostat - Serrated/Curved/14cm



Allgemein

|               |                            |
|---------------|----------------------------|
| Produktnummer | 13005-14                   |
| Kategorie     | Clamps & Clips > Hemostats |
| EAN           | 4251909503260              |



Finde 13005-14 auf [finescience.com](https://finescience.com)

Technische Informationen

|                     |                                 |              |      |
|---------------------|---------------------------------|--------------|------|
| Versandgewicht      | 30 g                            | Nettogewicht | 24 g |
| Paketgröße          | 1 Pc.                           | Maßeinheit   | EA   |
| Mindestbestellmenge | 1                               |              |      |
| Alloymaterial       | Surgical Stainless Steel (Inox) |              |      |
| Clampinglength      | 32 mm                           |              |      |
| Length              | 14 cm                           |              |      |
| Tipwidth            | 1.5 mm                          |              |      |
| Tips                | Serrated                        |              |      |
| Tipshape            | Curved                          |              |      |