



Fine Hemostat - Smooth/Straight/12cm



Allgemein

Produktnummer	13006-12
Kategorie	Clamps & Clips > Hemostats
EAN	4251909503277



Finde 13006-12 auf
finescience.com

Technische Informationen

Versandgewicht	25 g	Nettogewicht	18 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Clampinglength	19 mm		
Length	12 cm		
Tipwidth	1 mm		
Tips	Smooth		
Tipshape	Straight		