



Mosquito Hemostat - Serrated/Straight/12.5cm



Allgemein

Produktnummer	13008-12
Kategorie	Clamps & Clips > Hemostats Most Popular > Clamps & Clips
EAN	4061706012717



Finde 13008-12 auf
finescience.com

Technische Informationen

Versandgewicht	20 g	Nettogewicht	18 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Clampinglength	20 mm		
Length	12.5 cm		
Tipwidth	1.1 mm		
Tips	Serrated		
Tipshape	Straight		