



Mosquito Hemostat - Serrated/Curved/12.5cm



General Info

Product Number	13009-12
Category	Clamps & Clips > Hemostats Most Popular > Clamps & Clips
EAN	4061706012724



Find 13009-12 on
finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	18 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Clampinglength	20 mm		
Length	12.5 cm		
Tipwidth	1.1 mm		
Tips	Serrated		
Tipshape	Curved		