



Micro-Mosquito Hemostat - Serrated/Curved/12.5cm



General Info

Product Number	13011-12
Category	Clamps & Clips > Hemostats
EAN	4251909503307



Find 13011-12 on
finescience.com

Technical Info

Shipping Weight	25 g	Net Weight	21 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Clampinglength	20 mm		
Length	12.5 cm		
Tipwidth	1 mm		
Tips	Serrated		
Tipshape	Curved		