



Micro-Mosquito Hemostat - Serrated/Angled 90°/12.5cm



Allgemein

Produktnummer	13012-12
Kategorie	Clamps & Clips > Hemostats
EAN	4251909503314



Finde 13012-12 auf finescience.com

Technische Informationen

Versandgewicht	20 g	Nettogewicht	19 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Clampinglength	20 mm		
Length	12.5 cm		
Tipwidth	1 mm		
Tips	Serrated		
Tipshape	Angled 90°		