



Adson-Baby Hemostat - Serrated/Curved/14cm



Allgemein

Produktnummer	13013-14
Kategorie	Clamps & Clips > Hemostats
EAN	4251909503321



Finde 13013-14 auf
finescience.com

Technische Informationen

Versandgewicht	30 g	Nettogewicht	28 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Clampinglength	15.5 mm		
Length	14 cm		
Tipwidth	1.3 mm		
Tips	Serrated		
Tipshape	Curved		