



Adson-Baby Hemostat - Serrated/Curved/14cm



General Info

Product Number	13013-14
Category	Clamps & Clips > Hemostats
EAN	4251909503321



Find 13013-14 on
finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	28 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Clampinglength	15.5 mm		
Length	14 cm		
Tipwidth	1.3 mm		
Tips	Serrated		
Tipshape	Curved		