



Providence-Hospital Hemostat - Serrated/Straight/14.5cm



General Info

Product Number	13014-14
Category	Clamps & Clips > Hemostats
EAN	4251909503338



Find 13014-14 on
finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	27 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Clampinglength	15 mm		
Length	14.5 cm		
Tipwidth	1.5 mm		
Tips	Serrated		
Tipshape	Straight		