



Ultra Fine Hemostat - Smooth/Curved/12cm



These are the finest hemostats in our product line.

General Info

Product Number	13021-12
Category	Clamps & Clips > Hemostats
EAN	4251909503383



Find 13021-12 on
finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	20 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Clampinglength	19 mm		
Length	12 cm		
Tipwidth	0.6 mm		
Tips	Smooth		
Tipshape	Curved		