



Mosquito Hemostat - Serrated/Straight/12.5cm/Left-Handed



General Info

Product Number	13108-12
Category	Clamps & Clips > Hemostats Most Popular > Clamps & Clips
EAN	4251909503390



Find 13108-12 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	20 g	Net Weight	16 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Clampinglength	20 mm		
Length	12.5 cm		
Islefthanded	Yes		
Tipwidth	1 mm		
Tips	Serrated		
Tipshape	Straight		