



Mosquito Hemostat - Serrated/Curved/12.5cm/Left-Handed



General Info

|                |                                                             |
|----------------|-------------------------------------------------------------|
| Product Number | 13109-12                                                    |
| Category       | Clamps & Clips > Hemostats<br>Most Popular > Clamps & Clips |
| EAN            | 4251909503406                                               |



Find 13109-12 on  
[finescience.com](https://finescience.com)

Technical Info

|                        |                                 |                 |      |
|------------------------|---------------------------------|-----------------|------|
| Shipping Weight        | 20 g                            | Net Weight      | 15 g |
| Package Size           | 1 Pc.                           | Unit of Measure | EA   |
| Minimum Order Quantity | 1                               |                 |      |
| Alloymaterial          | Surgical Stainless Steel (Inox) |                 |      |
| Clampinglength         | 20 mm                           |                 |      |
| Length                 | 12.5 cm                         |                 |      |
| Islefthanded           | Yes                             |                 |      |
| Tipwidth               | 1.3 mm                          |                 |      |
| Tips                   | Serrated                        |                 |      |
| Tipshape               | Curved                          |                 |      |