



Standard Surgical Scissors - Straight/Blunt-Blunt/13cm/28mm



General Info

Product Number	14000-13
Category	Scissors & Micro Scissors > Surgical Scissors Most Popular > Scissors & Micro Scissors
EAN	4251909503413

Find 14000-13 on
finescience.com

Technical Info

Shipping Weight	40 g	Net Weight	35 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	42 mm		
Cuttingedge	28 mm		
Length	13 cm		
Serrated	No		
Tips	Blunt/Blunt		
Tipshape	Straight		