



Standard Surgical Scissors - Straight/Blunt-Blunt/14.5cm/30mm



Allgemein

Produktnummer	14000-14
Kategorie	Scissors & Micro Scissors > Surgical Scissors Most Popular > Scissors & Micro Scissors
EAN	4251909503420



Finde 14000-14 auf
finescience.com

Technische Informationen

Versandgewicht	40 g	Nettogewicht	40 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	45 mm		
Cuttingedge	30 mm		
Length	14.5 cm		
Serrated	No		
Tips	Blunt/Blunt		
Tipshape	Straight		