



Standard Surgical Scissors - Straight/Sharp-Blunt/14.5cm/30mm



General Info

Product Number	14001-14
Category	Scissors & Micro Scissors > Surgical Scissors Most Popular > Scissors & Micro Scissors
EAN	4251909503444



Find 14001-14 on
finescience.com

Technical Info

Shipping Weight	45 g	Net Weight	40 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	45 mm		
Cuttingedge	30 mm		
Length	14.5 cm		
Serrated	No		
Tips	Sharp/Blunt		
Tipshape	Straight		