



Standard Surgical Scissors - Straight/Sharp-Blunt/16.5cm/36mm



Allgemein

Produktnummer	14001-16
Kategorie	Scissors & Micro Scissors > Surgical Scissors Most Popular > Scissors & Micro Scissors
EAN	4061706016234



Finde 14001-16 auf finescience.com

Technische Informationen

Versandgewicht	65 g	Nettogewicht	59 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Bladelength	54 mm		
Cuttingedge	36 mm		
Length	16.5 cm		
Serrated	No		
Tips	Sharp/Blunt		
Tipshape	Straight		