



Standard Surgical Scissors - Straight/Sharp-Blunt/18.5cm/40mm



Allgemein

Produktnummer	14001-18
Kategorie	Scissors & Micro Scissors > Surgical Scissors Most Popular > Scissors & Micro Scissors
EAN	4061706016272



Finde 14001-18 auf finescience.com

Technische Informationen

Versandgewicht	75 g	Nettogewicht	73 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Bladelength	60 mm		
Cuttingedge	40 mm		
Length	18.5 cm		
Serrated	No		
Tips	Sharp/Blunt		
Tipshape	Straight		