



Standard Surgical Scissors - Straight/Sharp-Sharp/12cm/26mm



Allgemein

Produktnummer	14002-12
Kategorie	Scissors & Micro Scissors > Surgical Scissors Most Popular > Scissors & Micro Scissors
EAN	4061706016517



Finde 14002-12 auf finescience.com

Technische Informationen

Versandgewicht	30 g	Nettogewicht	27 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	39 mm		
Cuttingedge	26 mm		
Length	12 cm		
Serrated	No		
Tips	Sharp/Sharp		
Tipshape	Straight		