



Standard Surgical Scissors - Curved/Blunt-Blunt/12cm/26mm



Allgemein

Produktnummer	14003-12
Kategorie	Scissors & Micro Scissors > Surgical Scissors Most Popular > Scissors & Micro Scissors
EAN	4061706016524



Finde 14003-12 auf finescience.com

Technische Informationen

Versandgewicht	30 g	Nettogewicht	27 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	39 mm		
Cuttingedge	26 mm		
Length	12 cm		
Serrated	No		
Tips	Blunt/Blunt		
Tipshape	Curved		