



Standard Surgical Scissors - Curved/Blunt-Blunt/18.5cm/40mm



Allgemein

Produktnummer	14003-18
Kategorie	Scissors & Micro Scissors > Surgical Scissors Most Popular > Scissors & Micro Scissors
EAN	4061706016715



Finde 14003-18 auf finescience.com

Technische Informationen

Versandgewicht	70 g	Nettogewicht	68 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	60 mm		
Cuttingedge	40 mm		
Length	18.5 cm		
Serrated	No		
Tips	Blunt/Blunt		
Tipshape	Curved		