



Standard Surgical Scissors - Curved/Blunt-Blunt/18.5cm/40mm



General Info

Product Number	14003-18
Category	Scissors & Micro Scissors > Surgical Scissors Most Popular > Scissors & Micro Scissors
EAN	4061706016715

Find 14003-18 on finescience.com

Technical Info

Shipping Weight	70 g	Net Weight	68 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	60 mm		
Cuttingedge	40 mm		
Length	18.5 cm		
Serrated	No		
Tips	Blunt/Blunt		
Tipshape	Curved		