



Standard Surgical Scissors - Curved/Sharp-Blunt/14.5cm



Allgemein

Produktnummer	14004-14
Kategorie	Scissors & Micro Scissors > Surgical Scissors Most Popular > Scissors & Micro Scissors
EAN	4251909503505



Finde 14004-14 auf finescience.com

Technische Informationen

Versandgewicht	45 g	Nettogewicht	40 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Bladelength	45 mm		
Cuttingedge	30 mm		
Length	14.5 cm		
Serrated	No		
Tips	Sharp/Blunt		
Tipshape	Curved		