




Standard Surgical Scissors - Curved/Sharp-Blunt/14.5cm



General Info

Product Number	14004-14
Category	Scissors & Micro Scissors > Surgical Scissors Most Popular > Scissors & Micro Scissors
EAN	4251909503505



Find 14004-14 on
finescience.com

Technical Info

Shipping Weight	45 g	Net Weight	40 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Bladelength	45 mm		
Cuttingedge	30 mm		
Length	14.5 cm		
Serrated	No		
Tips	Sharp/Blunt		
Tipshape	Curved		