



Standard Surgical Scissors - Curved/Sharp-Sharp/12cm



General Info

Product Number	14005-12
Category	Scissors & Micro Scissors > Surgical Scissors Most Popular > Scissors & Micro Scissors
EAN	4061706016326

Find 14005-12 on
finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	23 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	39 mm		
Cuttingedge	26 mm		
Length	12 cm		
Serrated	No		
Tips	Sharp/Sharp		
Tipshape	Curved		