



Standard Surgical Scissors - Curved/Sharp-Sharp/13cm



General Info

Product Number	14005-13
Category	Scissors & Micro Scissors > Surgical Scissors Most Popular > Scissors & Micro Scissors
EAN	4251909503512

Find 14005-13 on
finescience.com

Technical Info

Shipping Weight	40 g	Net Weight	35 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	42 mm		
Cuttingedge	28 mm		
Length	13 cm		
Serrated	No		
Tips	Sharp/Sharp		
Tipshape	Curved		