



Standard Surgical Scissors - Serrated/Straight/Sharp-Blunt/14cm



Allgemein

Produktnummer	14007-14
Kategorie	Scissors & Micro Scissors > Surgical Scissors
EAN	4061706016760



Finde 14007-14 auf finescience.com

Technische Informationen

Versandgewicht	45 g	Nettogewicht	41 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	45 mm		
Cuttingedge	30 mm		
Length	14 cm		
Serrated	Yes		
Tips	Sharp/Blunt		
Tipshape	Straight		