



Standard Surgical Scissors - Serrated/Straight/Sharp-Blunt/14cm



General Info

Product Number	14007-14
Category	Scissors & Micro Scissors > Surgical Scissors
EAN	4061706016760



Find 14007-14 on finescience.com

Technical Info

Shipping Weight	45 g	Net Weight	41 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	45 mm		
Cuttingedge	30 mm		
Length	14 cm		
Serrated	Yes		
Tips	Sharp/Blunt		
Tipshape	Straight		