



Mayo Scissors - Curved/Blunt-Blunt/15cm



General Info

Product Number	14011-15
Category	Scissors & Micro Scissors
EAN	4251909503567



Find 14011-15 on
finescience.com

Technical Info

Shipping Weight	50 g	Net Weight	50 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	42 mm		
Cuttingedge	28 mm		
Length	15 cm		
Serrated	No		
Tips	Blunt/Blunt		
Tipshape	Curved		