



Mayo-Noble Scissors - Straight/Blunt-Blunt/17cm



General Info

Product Number	14014-17
Category	Scissors & Micro Scissors
EAN	4251909503604



Find 14014-17 on
finescience.com

Technical Info

Shipping Weight	90 g	Net Weight	86 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	45 mm		
Cuttingedge	30 mm		
Length	17 cm		
Serrated	No		
Tips	Blunt/Blunt		
Tipshape	Straight		