



Mayo-Noble Scissors - Curved/Blunt-Blunt/17cm



Allgemein

Produktnummer	14015-17
Kategorie	Scissors & Micro Scissors
EAN	4251909503611



Finde 14015-17 auf finescience.com

Technische Informationen

Versandgewicht	85 g	Nettogewicht	85 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	45 mm		
Cuttingedge	30 mm		
Length	17 cm		
Serrated	No		
Tips	Blunt/Blunt		
Tipshape	Curved		