



Metzenbaum Fino Scissors - Curved/Blunt-Blunt/13cm



General Info

Product Number	14019-13
Category	Scissors & Micro Scissors
EAN	4251909503673



Find 14019-13 on
finescience.com

Technical Info

Shipping Weight	25 g	Net Weight	19 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Bladelength	30 mm		
Cuttingedge	20 mm		
Length	13 cm		
Serrated	No		
Tips	Blunt/Blunt		
Tipshape	Curved		