



Metzenbaum Fino Scissors - Curved/Blunt-Blunt/14.5cm



General Info

Product Number	14019-14
Category	Scissors & Micro Scissors
EAN	4251909503680



Find 14019-14 on
finescience.com

Technical Info

Shipping Weight	25 g	Net Weight	24 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	33 mm		
Cuttingedge	22 mm		
Length	14.5 cm		
Serrated	No		
Tips	Blunt/Blunt		
Tipshape	Curved		