



Metzenbaum Fino Scissors - Curved/Blunt-Blunt/18cm



General Info

Product Number	14019-18
Category	Scissors & Micro Scissors
EAN	4251909503697



Find 14019-18 on
finescience.com

Technical Info

Shipping Weight	35 g	Net Weight	35 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	36 mm		
Cuttingedge	24 mm		
Length	18 cm		
Serrated	No		
Tips	Blunt/Blunt		
Tipshape	Curved		