



Iris Scissors - Curved/Sharp-Blunt/10cm



General Info

Product Number	14029-10
Category	Scissors & Micro Scissors > Fine Scissors
EAN	4251909503727



Find 14029-10 on
finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	16 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	25.5 mm		
Cuttingedge	17 mm		
Length	10 cm		
Serrated	No		
Tips	Sharp/Blunt		
Tipshape	Curved		