



Iris Scissors - Large Loops/Straight/Sharp-Sharp/10cm



Allgemein

Produktnummer	14040-10
Kategorie	Scissors & Micro Scissors > Fine Scissors
EAN	4251909503741



Finde 14040-10 auf finescience.com

Technische Informationen

Versandgewicht	20 g	Nettogewicht	18 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	24 mm		
Cuttingedge	16 mm		
Length	10 cm		
Feature	Large loops		
Haslargeloops	Yes		
Serrated	No		
Tips	Sharp/Sharp		
Tipshape	Straight		