



Standard Surgical Scissors - ToughCut®/Straight/Sharp-Blunt/13cm



General Info

Product Number	14054-13
Category	Scissors & Micro Scissors > Fine Scissors
EAN	4251909503765



Find 14054-13 on
finescience.com

Technical Info

Shipping Weight	40 g	Net Weight	36 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	42 mm		
Cuttingedge	28 mm		
Length	13 cm		
Feature	ToughCut®		
Serrated	Yes		
Tips	Sharp/Blunt		
Tipshape	Straight		