



Iris Scissors - ToughCut®/Straight/Sharp-Sharp/9cm



MOST
POPULAR



TOUGH
CUT



SERRATED

Allgemein

Produktnummer	14058-09
Kategorie	Scissors & Micro Scissors > Fine Scissors MedTech Line Most Popular > Scissors & Micro Scissors
EAN	4251909503772



Finde 14058-09 auf
finescience.com

Technische Informationen

Versandgewicht	15 g	Nettogewicht	12 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	24 mm		
Cuttingedge	16 mm		
Length	9 cm		
Feature	ToughCut®		
Serrated	Yes		
Tips	Sharp/Sharp		
Tipshape	Straight		