



Iris Scissors - ToughCut®/Straight/Sharp-Sharp/9cm



MOST
POPULAR



TOUGH
CUT



SERRATED

General Info

Product Number	14058-09
Category	Scissors & Micro Scissors > Fine Scissors MedTech Line Most Popular > Scissors & Micro Scissors
EAN	4251909503772



Find 14058-09 on
finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	12 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	24 mm		
Cuttingedge	16 mm		
Length	9 cm		
Feature	ToughCut®		
Serrated	Yes		
Tips	Sharp/Sharp		
Tipshape	Straight		