



Iris Scissors - ToughCut®/Straight/Sharp-Sharp/11.5cm



MOST  
POPULAR



TOUGH  
CUT



SERRATED

General Info

|                |                                                                                                       |
|----------------|-------------------------------------------------------------------------------------------------------|
| Product Number | 14058-11                                                                                              |
| Category       | Scissors & Micro Scissors > Fine Scissors<br>MedTech Line<br>Most Popular > Scissors & Micro Scissors |
| EAN            | 4251909503789                                                                                         |



Find 14058-11 on  
[finescience.com](https://finescience.com)

Technical Info

|                        |                                 |                 |      |
|------------------------|---------------------------------|-----------------|------|
| Shipping Weight        | 25 g                            | Net Weight      | 24 g |
| Package Size           | 1 Pc.                           | Unit of Measure | EA   |
| Minimum Order Quantity | 1                               |                 |      |
| Alloymaterial          | Surgical Stainless Steel (Inox) |                 |      |
| Bladelength            | 24 mm                           |                 |      |
| Cuttingedge            | 16 mm                           |                 |      |
| Length                 | 11.5 cm                         |                 |      |
| Feature                | ToughCut®                       |                 |      |
| Serrated               | Yes                             |                 |      |
| Tips                   | Sharp/Sharp                     |                 |      |
| Tipshape               | Straight                        |                 |      |