



Allgemein

Produktnummer	14060-10
Kategorie	Scissors & Micro Scissors > Fine Scissors Most Popular > Scissors & Micro Scissors
EAN	4251909503826



Finde 14060-10 auf
finescience.com

Technische Informationen

Versandgewicht	20 g	Nettogewicht	19 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	24 mm		
Cuttingedge	16 mm		
Length	10.5 cm		
Serrated	No		
Tips	Sharp/Sharp		
Tipshape	Straight		