

F · S · T

Fine Scissors - Straight/Sharp-Sharp/11.5cm



General Info

Product Number	14060-11
Category	Scissors & Micro Scissors > Fine Scissors
EAN	4251909503833



Find 14060-11 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	25 g	Net Weight	22 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Bladelength	24 mm		
Cuttingedge	16 mm		
Length	11.5 cm		
Serrated	No		
Tips	Sharp/Sharp		
Tipshape	Straight		