



Fine Scissors - Curved/Sharp-Sharp/11.5cm



General Info

Product Number	14061-11
Category	Scissors & Micro Scissors > Fine Scissors Most Popular > Scissors & Micro Scissors
EAN	4251909503864



Find 14061-11 on
finescience.com

Technical Info

Shipping Weight	25 g	Net Weight	21 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	24 mm		
Cuttingedge	16 mm		
Length	11.5 cm		
Serrated	No		
Tips	Sharp/Sharp		
Tipshape	Curved		