

F · S · T

Fine Scissors - Curved to Side/Sharp-Sharp/10.5cm



General Info

Product Number	14062-10
Category	Scissors & Micro Scissors > Fine Scissors
EAN	4251909503888



Find 14062-10 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	20 g	Net Weight	16 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloy material	Surgical Stainless Steel (Inox)		
Blade length	24 mm		
Cutting edge	16 mm		
Length	10.5 cm		
Serrated	No		
Tips	Sharp/Sharp		
Tip shape	Curved to side		