



Fine Scissors - Angled to Side/Sharp-Sharp/9cm



Allgemein

Produktnummer	14063-09
Kategorie	Scissors & Micro Scissors > Fine Scissors Most Popular > Scissors & Micro Scissors
EAN	4251909503895



Finde 14063-09 auf finescience.com

Technische Informationen

Versandgewicht	15 g	Nettogewicht	12 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	24 mm		
Cuttingedge	16 mm		
Length	9 cm		
Serrated	No		
Tips	Sharp/Sharp		
Tipshape	Angled to side		