



Stevens Scissors - Straight/Sharp-Sharp/11.5cm



Allgemein

Produktnummer	14064-11
Kategorie	Scissors & Micro Scissors
EAN	4251909503918



Finde 14064-11 auf finescience.com

Technische Informationen

Versandgewicht	25 g	Nettogewicht	21 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Bladelength	19.5 mm		
Cuttingedge	13 mm		
Length	11.5 cm		
Serrated	No		
Tips	Sharp/Sharp		
Tipshape	Straight		