



Stevens Scissors - Curved/Blunt-Blunt/11.5cm



Allgemein

Produktnummer	14067-11
Kategorie	Scissors & Micro Scissors
EAN	4251909503949



Finde 14067-11 auf
finescience.com

Technische Informationen

Versandgewicht	25 g	Nettogewicht	22 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	19.5 mm		
Cuttingedge	13 mm		
Length	11.5 cm		
Serrated	No		
Tips	Blunt/Blunt		
Tipshape	Curved		