



Stevens Scissors - Curved/Blunt-Blunt/11.5cm



General Info

Product Number	14067-11
Category	Scissors & Micro Scissors
EAN	4251909503949



Find 14067-11 on
finescience.com

Technical Info

Shipping Weight	25 g	Net Weight	22 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	19.5 mm		
Cuttingedge	13 mm		
Length	11.5 cm		
Serrated	No		
Tips	Blunt/Blunt		
Tipshape	Curved		