



Tissue Separating Scissors - Straight/Blunt-Blunt/11.5cm



Allgemein

Produktnummer	14072-10
Kategorie	Scissors & Micro Scissors > Fine Scissors
EAN	4251909503970



Finde 14072-10 auf finescience.com

Technische Informationen

Versandgewicht	25 g	Nettogewicht	19 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	24 mm		
Cuttingedge	16 mm		
Length	11.5 cm		
Serrated	No		
Tips	Blunt/Blunt		
Tipshape	Straight		