



Tissue Separating Scissors - Straight/Blunt-Blunt/11.5cm



General Info

Product Number	14072-10
Category	Scissors & Micro Scissors > Fine Scissors
EAN	4251909503970



Find 14072-10 on
finescience.com

Technical Info

Shipping Weight	25 g	Net Weight	19 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	24 mm		
Cuttingedge	16 mm		
Length	11.5 cm		
Serrated	No		
Tips	Blunt/Blunt		
Tipshape	Straight		