



Metzenbaum Mini Scissors - Curved/Blunt-Blunt/9cm



General Info

Product Number	14075-09
Category	Scissors & Micro Scissors
EAN	4251909504007



Find 14075-09 on
finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	12 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	18 mm		
Cuttingedge	12 mm		
Length	9 cm		
Serrated	No		
Tips	Blunt/Blunt		
Tipshape	Curved		